

ABT PVA CONCENTRATE



ABT PVA CONCENTRATE is a 10% concentrated solution used as a surfactant when added to Lactated Ringers or Hartmann Saline Solution. PVA's surfactant properties help prevent embryos from sticking to tubing, filters or dishes after collection. This product is an alternative to BSA concentrate but does not require refrigeration. It is packaged in a 20ml EVA bag and at 4ml/Liter, it is the same concentration of surfactant as ABT's Complete Flush.

This product was tested at Colorado State University during the 2018 breeding season where no embryos stuck to the petri dish when PVA Concentrate was used compared to 30% of embryos sticking to a petri dish where a surfactant was not present.

ABT products are all manufactured in a cleanroom environment and go thru strict quality control practices ensuring consistency and quality from lot to lot.

- Animal origin free & 100% traceable components
- Tight quality control standards equal consistent results
- Durable EVA packaging non-latex needle septum removable cap
- Made in the USA with two decades of experience
- latex needle septum removable cap
- Made in the USA with two decades of experience

<u>PRODUCT #</u>	<u>SIZE</u>
ABT020	20 ml

Storage: 2 - 30° C
Shelf Life: 18 months from date of manufacture

ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H
PULLMAN, WA 99163 USA
Phone: 509-592-8144
www.abt360llc.com
alyssa@abt360llc.com
facebook.com/ABT360
instagram.com/ABT_360