

ABT FREEZE (10% GLYCEROL)



ABT FREEZE (GLYCEROL) is used for freezing embryos for storage and contains 10% Glycerol and 0.1M Sucrose to help cells from over expanding. Rehydration is required when thawing.

ABT's Freeze starts with ultra-purified Type 1 water. All ABT products are manufactured in a cleanroom environment and go thru strict quality control practices, ensuring consistency and quality from lot to lot.

- 10% Glycerol and 0.1M Sucrose
- Includes Gentamycin as an antibiotic
- 100% traceable components including certified New Zealand sourced BSA
- Tight quality control standards equal consistent results
- Available in durable EVA bags and 7ml LDPE Single Use Dropper Bottle
- Made in the USA with two decades of experience

<u>PRODUCT #</u>	<u>SIZE</u>
ABT032	50 mL Bag
ABT232	20 mL Bag
ABT732	7 mL Bottles, Pack of 5

Storage: 2 - 8° C
Shelf Life: 18 months from date of manufacture

ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H
PULLMAN, WA 99163 USA
Phone: 509-592-8144
www.abt360llc.com
alyssa@abt360llc.com
facebook.com/ABT360
instagram.com/ABT_360